# STORZ

## KARL STORZ — ENDOSKOPE

203020 20-1

EQUIMAT®

203020 20-1

EQUIMAT®

203020 20-1

EQUIMAT®

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Instruction Manual

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**Important information  
for users of KARL STORZ  
instruments**

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.




This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the EQUIMAT®. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

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## Controls, displays, connectors, and their uses

- ① Power switch
- ② M button
- ③ ± buttons for setting the alarm limit values
- ④ Operating mode symbol:  
volume difference measurement
- ⑤ Operating mode symbol:  
volume difference flow measurement
- ⑥ Display
- ⑦ Reset button
- ⑧ Mute button for acoustic alarm
- ⑨ Label strip (replaceable)
- ⑩ Socket for weighing cell
- ⑪ Potential equalization connector
- ⑫ Power fuse holder
- ⑬ Power socket
- ⑭ SCB connectors \*
- ⑮ Weighing cell
- ⑯ Suspension holder, consisting of:
  - ⑰ Irrigation-liquid bag holder
  - ⑱ Square tube
  - ⑲ Tube holder
  - ⑳ Bottle holder

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\* The KARL STORZ-SCB® interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

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
	Read the instructions carefully before operating the equipment!
	Connection for weighing cell
	Reset
	Accoustic alarm off
	Volume difference measurement
	Volume difference flow measurement
	Potential equalization connector
	Alternating current

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	<b>DANGER:</b> Risk of explosion if used in the presence of flammable anesthetics.
	<b>CAUTION:</b> To avoid the risk of electrical shock, do not remove cover. Refer servicing to qualified service personnel.  Keep out of reach of patients. Do not store liquids on or above the unit.

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Photographs of the unit .....	A 1
Controls, displays, connectors, and their uses .....	A 2
Symbols employed .....	A 3

## **General information**

Device description .....	2
--------------------------	---

## **Safety instructions**

Warnings and cautions .....	3
Intended use .....	6
User qualification .....	6
Safety precautions at the site of installation .....	7
Safety precautions when using the unit .....	7
Safety devices .....	8

## **Installation and operating instructions**

Unpacking the equipment .....	9
Basic equipment .....	9
Installing and connecting up the unit .....	9
Operating the unit .....	16
Volume difference display mode .....	17
Volume difference flow display mode .....	17
Actual value/alarm limit display .....	17
Alarm signals .....	18
Changing the dispenser bag .....	19

## **Maintenance**

Fuse replacement .....	20
Cleaning and sterilization .....	21
Maintenance .....	26
Servicing and repair .....	26
Repair program .....	27
Limitation of liability .....	28
Manufacturer's warranty .....	28

## **Technical description**

Troubleshooting .....	29
Technical data .....	30
Technical documentation .....	31
General circuit diagram .....	32

## **Spare parts, accessories**

Spare parts list .....	33
------------------------	----

## **Appendix**

Cleaning agents and disinfectants .....	34
Electromagnetic Compatibility (EMC) Information for EQUIMAT® .....	36

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## Device description

To minimise the risk of excessive intravasation of the irrigation fluid, the EQUIMAT® permits monitoring of the volume difference between inflowing and outflowing irrigation liquid.

For this purpose, the EQUIMAT® has 2 operating modes:

- volume difference display, which indicates the difference between inflowing and outflowing quantity of irrigation liquid, and
- volume difference flow display, which indicates the difference between inflowing and outflowing irrigation liquid per unit of time (ml/min).

Alarm limit values can be defined for both operating modes and when these are reached an alarm signal is given to warn the user.

The EQUIMAT® is designed for use in conjunction with the KARL STORZ suction/irrigation pump range, but can also be used with corresponding instruments from other manufacturers.

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## Warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words, **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

**Warning:** A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

**Caution:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

**Note:** A Note indicates special information about operating the product, or clarifies important information.

**Warning:** Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

**Warning:** The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

**Warning:** Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cp. IEC 60601-1-1).

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**Warning:** Test this equipment prior to each surgical procedure to ensure that it functions correctly.

**Warning:** Risk of explosion if used in the presence of flammable anesthetics (see "Safety precautions at the site of installation", page 7).

**Warning:** To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

**Warning:** Grounding reliability can only be achieved when the equipment is connected to "Hospital Only" or "Hospital Grade" outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

**Warning:** The electrical installation of the operating room must comply with the applicable IEC requirements.

**Warning:** Keep out of reach of patients.

**Warning:** Always unplug the unit before doing any work on it (e. g. cleaning).

**Warning:** The EQUIMAT® has been programmed at the factory with physiologically sensible alarm limit values. If these are exceeded, the instrument gives an acoustic and visual warning. These alarm limit values can be changed by the user. It is the physician's responsibility to adjust these values according to the particular application and patient and to interpret the medical significance of the values measured.

**Warning:** The indicated volume difference between inflowing and outflowing irrigation liquid can be falsified, for example by leaks in the tubing system or by losses of liquid when an instrument is changed. In this case the values measured by the EQUIMAT® give the impression of a higher intravasation than is actually the case.

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**Caution:** Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

**Caution:** Before connecting the unit to the electrical supply lines, verify that the line voltage stated on the unit's identification plate agrees with that of the electrical supply lines to be used.

**Caution:** Use only fuses of the correct rating.

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**Note:** This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

**Note:** Make certain to observe national regulations/laws regarding disposal.

**Note:** Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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### **Intended use**

The EQUIMAT® is designed for use in conjunction with KARL STORZ suction/irrigation pumps and also suitable units from other manufacturers. It measures the volume balance between inflowing and outflowing irrigation liquid during endoscopic interventions.

Use of the unit in fields other than those indicated above is not allowed for safety reasons.

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The EQUIMAT® may only be used with accessories which have been designated by KARL STORZ as suitable for the unit or the safe use of which is proven.

Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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
### **User qualification**

The EQUIMAT® may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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### **Safety precautions at the site of installation**

The unit may only be used in medical rooms of which the installations conform to applicable national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the demarcated hazard zone. Examples of such substances are: anesthetic ether (diethyl ether, cyclopropane) as well as combustible, volatile skin cleansers and skin disinfectants which may create an explosive atmosphere (e.g. detergent ether, petroleum ether).

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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### **Safety precautions when operating the unit**

It is the user's responsibility to make sure the unit is safe and operates properly before using it. During treatment using the EQUIMAT® the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

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## Safety features

### Self-test procedure

Immediately after switching on, the EQUIMAT® carries out a self-test, in which the unit and the weighing cell are checked. If the result of the self-test is negative, the keypad is blocked and an error message appears in the display. Correct functioning of the weighing cell is also checked regularly during operation.

### Zero point balancing

To obtain correct measurement, the unit must be first set to zero by pressing the RESET button. If this reset is not carried out, the keypad remains blocked.

### Alarm limit values

The unit is programmed at the factory with 2 preset alarm limit values. When these are reached, an alarm is triggered (intermittent acoustic signal and flashing display).

The alarm is triggered at the following limit values:

- volume difference measurement: 500 ml
- volume difference flow measurement: 100 ml/min

The alarm limit value can be redefined by the user at his own discretion.

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### Warning of positive liquid balance

If more liquid is drained than is flowing into the patient, this is an indication that significant quantities of body fluid are being removed. In this case a warning signal appears in the form of a flashing plus sign in the display.

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### Unpacking the equipment

Carefully unpack the EQUIMAT® and remove the unit and its accessories from their packaging.  
Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately.  
If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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### Basic equipment

- 1 EQUIMAT® 203020 20-1
- 1 Weighing cell with connecting cord, 383321 30
- 1 Suspension holder 203020 31
- 1 SCB connecting cord 200901 70, length 100 cm
- 1 Power cord 400 A
- 1 Instruction manual 96126001 d
- 1 Abridged instructions (self-adhesive)

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
### Installing and connecting up the unit

**Note:** The EQUIMAT® including the accessories connected may be used only in medical facilities having electrical installations conforming to the applicable national, state, and local electrical codes.

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**Warning:** This unit is not intended for use in hazardous zones. Do not operate the unit within demarcated hazard zones while explosive anesthetic gases are in use.

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**Fitting a label strip to the front panel**

Undo the left-hand front side cover with a screwdriver.

**Note:** The appropriate controls/indicators on the unit are designated on the label strip.

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Remove any unmarked label strip. Select a label strip ⑨ (currently 6 languages available) and slide it between the plastic and the metal front. Fold the strip along the broken line to an angle of 90° and approx. 2 mm further on (toward the end of the strip) fold it back to an angle of 90° so that the side cover can be refitted without any obstructions. Replace the side cover.

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Install the unit in the immediate proximity of the suction/irrigation pump, preferably on the pump itself.

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The unit is equipped with a connector ⑩ for attaching a ground line. The unit's ground line should be installed by a qualified electrician.

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## Installation and operating instructions

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Before plugging in the unit, make sure that the voltage on the nameplate corresponds to the voltage of the local powerline.

Connect power cord. Insert power plug into power socket 13 as far as it will go.

**Warning:** Insert the power cord only into electrical outlets located outside areas subject to explosion hazards.

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### KARL STORZ-SCB®

**Note:** To avoid mistakenly pulling out the SCB connecting cable, the SCB connector possesses a protection device.

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets 14. Connect the other end of the cable to the KARL STORZ-SCB® (KARL STORZ Communication Bus) control unit or other SCB units (see KARL STORZ-SCB® System Instruction Manual).

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
### Assembling the suspension holder

Unscrew and remove the hex socket bolts from the bottle holder 20 using the tool enclosed.

Insert the bottle holder 20 into the square tube at the end facing the triangle that is not filled in. (The 5 support elements are located at the same end as the printed arrow markings. Refer to diagram of suspension holder on page A1).

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**Note:** The correct end of the square tube has a notch into which a pin on the bottle holder 20 fits, so that the parts only fit together correctly in one position.

Secure bottle holder to square tube by screwing the bolts tight.

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
Fit tube holder (19) onto square tube (18) (the filled-in black arrow markings printed on the tube holder and on the square tube must be aligned. Refer also to picture of suspension holder on page A1).

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Unscrew and remove the hex socket bolts from the irrigation-liquid bag holder (17) with the tool enclosed. Insert the irrigation-liquid bag holder into the square tube (18) (the hooks of the bag holder (17) are located at the same end as the printed arrow markings. Refer to picture of suspension holder on page A1).

Secure the irrigation-liquid bag holder to the square tube and screw the bolts tight.

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Screw the weighing cell (15) tightly onto the stand holder for irrigation liquid. Make sure that the weighing cell is horizontal.

**Note:** When selecting the stand system it must be ensured that the expected maximum load of 30 kg system weight can be suspended without any risk. KARL STORZ recommends the mobile stand 20 0200 24.

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Plug the connecting cord of the weighing cell into the socket (10) of the EQUIMAT®.

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
Attach the suspension holder ⑯ to the weighing cell ⑮.

**Note:** The weighing cell must be connected and the suspension holder must be fitted before switching on the unit.

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Screw the tube holder ⑲ to the irrigator rod so that the weighing cell and suspension holder are aligned parallel to the irrigator rod (not at an angle).

**! Caution:** Risk of incorrect measurements!  
Make sure that the holder is located within the two limit arrows printed on the square tube.

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Attach the bottle holder to the appropriate retainer of the suspension holder ⑯. Fit the cap onto the suction bottle and insert the bottle into the holder.

**Note:** Make sure that the cap is correctly fitted and tightly closed, as otherwise sufficient suction force cannot be achieved.

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### Using more than one suction bottle – Mobile Stand

**Note:** When using the Mobile Stand, the bottles must be fitted symmetrically at the back and front (to ensure stability).

With STORZ bottles: one at the front, two at the back

With plastic bottles: three at the front, two at the back

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**Using more than one suction bottle  
– Video Trolley**

**Note:** When using a video trolley, the suction bottles must be fitted symmetrically on the front of the suspension holder.

With STORZ bottles: up to three

With plastic bottles: up to five

If 5 plastic bottles are used on the video trolley, the optional inserts 20 3020 32 must be fitted.

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**Note:** If suction bottles are arranged in series, only the “last” bottle may be fitted with an over-flow protection device. The float must be removed from the other bottles.

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
Two or more suction bottles must be joined together with a short connecting tube (art. no. 20 3002 42).

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**Preparing the suction pump**

Fit the suction tube (more than one where necessary).

Connect the suction tube to the suction pump (Endomat®) by means of LUER lock connection (male).

**Note:** To protect the suction pump, we recommend using a hydrophobic filter. If necessary, refer to the instruction manual for the suction pump.

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Secure the suction tube (more than one where necessary) to the tube holder at the bottom (and top).

**Note:** The suction tube must have a decoupling loop between the suction bottle and the tube holder.

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Never secure the tubes like this!

**Warning:** Direct fitting of the tubes (without a loop) leads to measuring errors and must be prevented under all circumstances!

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
Attach irrigation-liquid bag symmetrically.

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### Prepare the irrigation pump

Fit the irrigation tube as described in the instruction manual for the irrigation pump.

Secure irrigation tubing to the tube holder at the top.

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## Operating the unit

Switch on the unit at the power switch ①. The unit now carries out a self-test. If no error is found, the weight of the suspension holder and the fitted irrigation-liquid bags and suction bottles appears on the display ④.

**Note:** Make sure that the maximum permissible weight of 30 kg is not exceeded. For weights larger than 20 kg (> 19,999 g) the displayed value changes. The units digit disappears and values in the 1,000s are displayed with a decimal point. For example, 20,000 ml is shown as 20.00 ml. The largest possible value that can be displayed is  $\pm 32,760$ , which means  $\pm 32.76$  on the display.

If an error is detected, the letter E and an error code appear in the display. In this case the keypad remains blocked until the error has been rectified.

**Note:** To reduce measuring errors, we recommend removing the air from the tube system before the zero adjustment (Reset button).

Press the reset button ⑦ to set the volume display to 0.

**Note:** The keypad remains blocked until the reset has been carried out.

Then check calibration of the weighing cell. For this, also attach a known weight (e.g. 3,000 ml irrigation-liquid bag) to the weighing cell or remove a bag which has already been attached. Check the change in volume display. It must be equal to the test weight ( $\pm 30$  ml) and have the appropriate sign + or -. The unit is now ready for operation.

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**Volume difference display  
operating mode**

Every time the EQUIMAT® is switched on it initializes in volume difference display mode (the corresponding operating mode symbol ④ glows).

In this operating mode the EQUIMAT® shows the volume difference between the quantity of irrigation liquid flowing into the body and the quantity being removed again.

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**Volume difference flow display  
operating mode**

The volume difference flow display indicates the irrigation liquid loss as a function of time, i.e. the difference between the flow of irrigation liquid entering the body and the drainage flow in ml/min. To activate the volume difference flow display, the M button ② must be pressed briefly (< 2 sec). The corresponding operating mode symbol ⑤ lights up.

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**Actual value/alarm limit display**

The EQUIMAT® is programmed at the factory with an alarm limit value in both volume difference display mode and volume difference flow display mode. When the limit value is reached an alarm is triggered. The limit value can be displayed in any operating mode by pressing the M button ② for longer than 2 seconds. The alarm limit value display is identified by a flashing decimal point on the right.

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The operator can change the alarm limit value at his own discretion using the  $\pm$  buttons ③. If the M button ② is pressed again briefly, the actual values are displayed again. The display also returns to the actual values automatically 10 seconds after the M button or one of the  $\pm$  buttons was last pressed.

**Note:** Changes to the alarm limit values are not stored. The unit initializes with the values preset at the factory every time it is switched on.

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## Alarm signals

If the alarm limit value is reached in the active operating mode, an alarm is given.

**Note:** The user must explicitly select the operating mode which he wishes to use.

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## Alarm function in the volume difference mode

If the alarm limit value is exceeded, an acoustic alarm is triggered after a latency of approx. 10 sec. (to suppress measurement artefacts), and the display begins to flash.

The acoustic alarm can be acknowledged (deactivated) by pressing the mute button ⑧. The display flashes for as long as the limit value is exceeded.

Five minutes after the acoustic alarm is muted, it sounds again (reminder) even if the volume difference has not increased any further.

If the new alarm limit (+125 ml) is exceeded again, the volume difference alarm sounds immediately.

This process can be repeated any number of times.

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### Alarm function in the volume difference flow mode

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The acoustic alarm can be acknowledged (deactivated) by pressing the mute button ⑧.

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If a volume increase is measured instead of the expected volume loss (e.g. because significant quantities of blood are drained together with the irrigation liquid), the sign in the display ④ changes from – to + and starts to flash.

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### Changing the dispenser bag


To change the dispenser bag during measurement, press the reset button ⑦ briefly (< 3 sec). The value „0000“ flashes in the display ⑥, but the current volume difference remains stored. After the irrigation liquid bag has been changed, measuring is resumed by again briefly pressing the reset button ⑦.

**Note:** If it takes longer than 30 sec. to change the dispenser bag, an alarm signal sounds which can be turned off by pressing the reset button after the bag has been changed.

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### Resetting the volume display during measurement

To set the volume display to zero during operation, the reset button ⑦ must be held down for longer than 3 seconds. Immediately after pressing the reset button, the value „0000“ flashes in the display and after 1 second an acoustic warning signal is given. After 3 seconds the display is set to zero and the previous measurement is erased.

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### Fuse replacement

Switch off the unit and remove the power plug from the electrical outlet.

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Remove the power fuse holder ⑫ with a screwdriver or other tool.

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**! Caution:** Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

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Replace and the power fuse holder ⑫.  
Reconnect the power cord.  
Test the unit for proper operation.

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## Cleaning and sterilization

**Note:** The instructions for “Cleaning, Sterilization and Care of KARL STORZ Instruments” must be followed. There the procedures for Cleaning, Disinfection and Sterilization are explained in detail.

**Warning:** Always pull out power plug before cleaning!

**Caution:** Avoid allowing fluids to enter the unit.

The exterior of the unit and the weighing cell can be cleaned with a conventional surface disinfectant.

**Warning:** Wear protective gloves, clothing and a face mask for cleaning of contaminated tubing.

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## Cleaning and sterilization of tubing

**Warning:** The tubing should be thoroughly cleaned and sterilized according to validated infection control procedures prior to use and subsequent reuse.

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**Warning:** High-level disinfection is NOT recommended for any tubing to be used for laparoscopy or gynecology. In such cases sterilization is recommended. The disinfection referred to here is not the disinfection which is performed for the protection of personnel.

**Warning:** Discard DISPOSABLE tubing, locks and bacterial hydrophobic filters.

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### **Preparation for cleaning and sterilization**

Place the tubing set in a container filled with a cleaning solution (per manufacturer's instructions) immediately after use to prevent contaminants from coagulating and drying onto the tubing set.

### **Water Quality Requirements**

Demineralized water is recommended for cleaning and rinsing and preparing the cleaning/disinfectant solution.

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## **Cleaning**

### **Manual cleaning**

Thoroughly clean the bottles and the overflow protection valve. Clean the channels with a brush and then rinse out for one minute. The rinsing water must no longer contain any dirt particles. Ensure that the floats of the overflow protection unit can move freely and do not stick due to residues.

The bottle caps, including the overflow protection valve can be steam sterilized or gas sterilized or be immersed in solution.

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## **Chemical Disinfection**

Use high level disinfectant solutions which are specially approved for endoscopes from KARL STORZ (see Appendix, pages 34, 35).

**Caution:** Follow the solution manufacturer's instructions paying close attention to proper dilution and recommended minimum exposure times. Pro-longed exposure may result in corrosion or material damage.

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**Cleaning the tubing set manually**

**Warning:** Wear protective gloves, clothing and a face mask when cleaning a contaminated tubing set.

1. Thoroughly rinse the tubing set to remove all coarse dirt particles.
2. Completely immerse the tubing set (be sure that the inside of the tubing set is completely filled with cleaning solution) in a cleaning solution (per manufacturer's instructions) prepared with demineralized water. KARL STORZ advises you not to use any detergents other than those approved (see Appendix).
3. With the help of a syringe, rinse out the tubes three times. Remove contaminants with brushes, sponges, soft cloths or a cotton applicator. *Cleaning accessories are available from KARL STORZ (see catalog).*
4. Cleaning brushes must be cleaned and disinfected immediately.
5. After disinfection, triple-rinse all tubing sets with demineralized water, for a minimum of one minute for each rinse. The rinse water should be discarded at the end of each rinse, as it will be contaminated with the cleaning solution. Rinse the tubing set thoroughly outside and inside in order to remove any residue which could interfere with sterilization. A cleaning pistol (Art. no. 27660) with a fine jet attachment is useful for rinsing the tubing sets.
6. Dry the tubing set with a lint-free soft cloth or sterile compressed air. The inner lumen in particular must be absolutely dry. A cleaning pistol (Art. no. 27660) with a fine jet attachment is useful for drying the tubing sets with compressed air.
7. After cleaning, inspect the tubing for cleanliness and damage.

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**Caution:** Follow the solution manufacturer's instructions paying close attention to proper dilution and recommended minimum exposure times. Prolonged exposure may result in corrosion or material damage.

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### Cleaning by machine

To clean by machine, connect the tubing to the retaining mechanisms in such a way that the instrument channels can be completely rinsed out.

The tubing can be cleaned and disinfected chemically or thermally.

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**Note:** It is recommended that the entire tubing be checked for leaks after cleaning. As a result it is possible to detect a leak early, through which liquid could escape and enter the unit.

To check for leaks use leak tester 13242 XL for example. For this purpose connect one end to the leak tester and block off the other end, with your finger for example.

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### Sterilization

Routine steam sterilization at 134°C (+3°C) is recommended before initial and subsequent use of all tubing.

Clean and dry the tubing as described in the section 'Cleaning'. Then place the tubing in suitable packaging.

**Warning:** Sterilization is only possible on clean surfaces.

**Warning:** The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and validated.

**Warning:** Any deviations from the recommended parameters for sterilization should be validated by the user.

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### Steam Sterilization

KARL STORZ has validated the following sterilization method:

Fractionated prevacuum (steam sterilization)

**Warning:** Tubing sets may only be sterilized using the fractionated prevacuum procedure.

Trays should be positioned in the sterilizer so that there is adequate circulation and penetration of steam, air removal and condensate drainage.

Load the sterilizer according to its instruction manual.

At the completion of the steam sterilization cycle, all instruments should remain untouched until adequately cooled.

**Warning:** Burn hazard! Parts are hot after steam sterilization. Allow them to cool down.

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### Fractionated prevacuum procedure

The fractionated prevacuum procedure consists of four phases.

- **Conditioning phase**  
The conditioning phase removes air from the chamber by pulling a vacuum up to four times. Steam is then blown in.
- **Sterilization phase**  
Sterilization takes place at a temperature of 134°C (+3°C) and an Operating pressure ( $p_b$ ) of 200 kPa (2 bar/29 psi) over a period of 5 minutes (minimum time).
- **Evacuation phase**  
The steam is released.
- **Drying phase**  
Drying takes place under renewed application of a vacuum over a period of approx. 5–20 minutes.

Validated parameters for the fractionated prevacuum method:

Temperature:	134°C (+3°C)
Pressure ( $p_b$ ):	200 kPa (2 bar/29 psi)
Exposure time:	5 minutes for all instruments.

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
## Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

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## Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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## Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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**Repair program**

For fiberscopes and equipment, individual repair is necessary. Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Dr. Karl-Storz-Str. 34  
78532 Tuttlingen.

In other countries please contact your local KARL STORZ branch or authorized dealer.

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**Important information**

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/equipment to the sender. Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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## Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

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## Manufacturer's warranty

For a period of **two years** after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply. Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG  
Postfach 230  
78503 Tuttlingen/Germany

**Note:** The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Dr. Karl-Storz-Str. 34  
78532 Tuttlingen/Germany

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Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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### Troubleshooting

**Warning:** Always unplug the unit beforehand!

*Symptom:*  
Complete failure of the unit.  
Display does not light up.

*Possible causes:*  
No power from the power line.  
Defective power fuse.

*Remedy:*  
Check that there is electricity to the wall outlet.  
Change fuses as described in the instruction manual.  
Make sure you fit the correct type of fuse.

*Symptom:*  
Error code E.

*Possible causes:*  
Volume difference flow measuring range of the weighing cell has been exceeded (>20 l/min.).

*Remedy:*  
Reduce volume difference flow.

*Symptom:*  
Error code E 1.

*Possible causes:*  
Weighing cell defective or not connected.

*Remedy:*  
Check connection of weighing cell or fit new weighing cell.

All error codes > E 1 entail an internal fault and the instrument must be inspected by an authorised repair service.

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**Standard compliance (for 203020 20-1)**

**According to IEC 60601-1, UL 2601.1, CAN/CSA C22.2 No. 601.1-M90:**

Type of protection against electric shocks:  
Protection Class I

**According to IEC 60601-1-2:**

Please read the Electromagnetic Compatibility Information in the appendix (p. 36-48).

**Directive compliance (for 203020 20-1)**

**According to Medical Device Directive (MDD):**

This medical product belongs to Class I

This medical product bears the CE mark in accordance with MDD 93/42/EEC. A code number after the CE mark indicates the responsible notifying body.

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**Согласно Medical Device Directive (MDD):**

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### Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated as repairable by their respective manufacturers.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

*We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.*

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**Type**

- 1 = Cleaner  
2 = Disinfectant  
3 = Suitable for ultrasound baths  
A = Not suitable for flexible endoscopes  
(according to manufacturers' instructions)

**Note:** You will also find the current list of approved products on our Web site ([www.karlstorz.com](http://www.karlstorz.com)).

**Note:** When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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**Chemicals for the preparation of telescopes and endoscopic instrument sets**

**Manual cleaning and disinfection/ultrasound**

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, prismatic light deflectors. For motors and their handles the appropriate manual provides detailed information.

The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to proper dilution and recommended immersion times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility (Note: For video endoscopes please refer to our separate list of approved chemicals.):

Ecolab GmbH & Co. OHG	Sekusept forte	2, 3
	Sekusept Plus	1, 2, 3
	Sekusept Pulver Classic	1, 2, 3
	Sekusept Pulver Classic + Aktivator	1, 2, 3
Esteer GmbH	Ultrasemit Konz.	2, 3
Fresenius AG	Afid	1, 2
	Afid plus (neu)	1, 2
	Sporcid FF	1, 2
Indeba Industria E Commercio LTDA	Letahdeido	2
Jose Collado S.A.	Darodor 4000 Liquido	1
	Darodor 9000	2
	Darodor Sinaldehyd 2000	2
Lysoform Dr. Rosemann GmbH	Aldasan 2000	1, 2
	Almyrol	1, 2
	Desoform	1, 2
	Lysoformin 3000	1, 2
Merz+Co.	Edisonite Super	3
	Mucadont-IS	2, 3
	Mucadont-Zymaktiv	1, 3
	Mucocit-T	1, 2, 3
Orochemie	A 10	2
	ID 210	2
Promagent AB	Wavacide	2
Schülke & Mayr GmbH	Gigasept	2, 3
	Gigasept FF	2, 3
	Lysetol AF	1, 2, 3
	Lysetol V	2
Schuhmacher, Dr. GmbH	Descoton forte	2
	Perfektan Endo	2
Steris	EnzyCare 2	1, 3
	neodisher LM 2	1
Weigert, Dr. GmbH & Co.	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher medizym	1, 3
Whiteley Industries PTY.Ltd	Aidal Plus	2

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**\* Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set.

For **flexible endoscopes**, only immerse the examination sheath, **not** the housing/handle part.

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**Type**

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- 4 = Neutralizer
- 5 = Rinse aid

**Note:** You will also find the current list of approved products on our Web site ([www.karlstorz.com](http://www.karlstorz.com)).

**Note:** When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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**Machine cleaning and disinfection**

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack. For motors and their handles the appropriate manual provides detailed information.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Schuhmacher, Dr. GmbH	Thermoton Endo	2
Steris	Instru-Klenz	1, 3
	neodisher FA	1
Weigert, Dr. GmbH & Co.	neodisher FA forte	1, 3
	neodisher mediclean	1
	neodisher mediclean forte	1
	neodisher mediklar	5
	neodisher medizym	1
	neodisher Septo DA	2
	neodisher Septo DN	2

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## Appendix

### Electromagnetic Compatibility (EMC) Information for EQUIMAT®

**WARNING:** Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The EQUIMAT® model 203020 20-1 corresponds to EN/IEC 60601-1-2 :2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC. These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The EQUIMAT® model 203020 20-1 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

**NOTE:** The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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## Appendix

### Electromagnetic Compatibility (EMC) Information for EQUIMAT®

**WARNING:** The EQUIMAT® model 203020 20-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

**WARNING:** The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

**WARNING:** The use of accessories, transducers and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the EQUIMAT® model 203020 20-1. The accessories, transducers and cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories, transducers and cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Firmenstempel/Unterschrift:

# WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

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To be completed by the supplier/importer:

Company stamp/signature:

To be filled out by instrument owner:

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KARL STORZ Endoscopy China Ltd.  
Hong Kong Representative Office  
Unit 1601, Chinachem Exchange Square  
1 Hoi Wan Street, Quarry Bay, Hong Kong,  
People's Republic of China  
Phone: +8 52 28 65 2411  
Fax: +8 52 28 65 4114  
E-Mail: [inquiry@karlstorz.com.hk](mailto:inquiry@karlstorz.com.hk)

KARL STORZ Endoscopy China Ltd.  
Beijing Representative Office  
Room 610, China Life Tower  
No. 6, Chaowai Street  
Beijing, 100020, People's Republic of China  
Phone: +86 10 8525 3725  
Fax: +86 10 8525 3728  
E-Mail: [inquiry@karlstorz.com.hk](mailto:inquiry@karlstorz.com.hk)

KARL STORZ Endoscopy China Ltd.  
Unit 3901-3904, Tower 1 Grand Gateway  
No. 1 Hong Qiao Road  
Shanghai, 200030, People's Republic of China  
Phone: +86 21 6113-1188  
Fax: +86 21 6113-1199  
E-Mail: [inquiry@karlstorz.com.hk](mailto:inquiry@karlstorz.com.hk)

KARL STORZ Endoscopy China Ltd.  
Chengdu Representative Office  
F-5, 24/F, Chuanxing Mansion  
No. 18 Renming Road South  
Chengdu, Sichuan, 610016, People's Rep. of China  
Phone: +86 28 8620-0175  
Fax: +86 28 8620-0177  
E-Mail: [inquiry@karlstorz.com.hk](mailto:inquiry@karlstorz.com.hk)

KARL STORZ Endoscopy China Ltd.  
Guangzhou Representative Office  
Room 1119-20, Dongshan Plaza  
69 Xianlie Road Middle, Dongshan District,  
Guangzhou, Guangdong, 510095,  
People's Rep. of China  
Phone: +86 20 8732-1281  
Fax: +86 20 8732-1286  
E-Mail: [inquiry@karlstorz.com.hk](mailto:inquiry@karlstorz.com.hk)

KARL STORZ Endoscopy Japan K. K.  
Bois Hongo Building 6Fl, 3-42-5 Hongo  
Bunkyo-ku, Tokyo 113-0033, Japan  
Phone: +81 3 5802-3966  
Fax: +81 3 5802-3988  
E-Mail: [info@karlstorz.co.jp](mailto:info@karlstorz.co.jp)

## **KARL STORZ GmbH & Co. KG**

Mittelstraße 8  
78532 Tuttlingen

Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 (0)7461 708-0  
Telefax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.de](mailto:info@karlstorz.de)  
Web: [www.karlstorz.com](http://www.karlstorz.com)